

Needle/Glass Trap

CLEANROOM VACUUM CLEANER ATTACHMENT
FOR NEEDLE AND GLASS RECOVERY
SUITABLE FOR PHARMACEUTICAL APPLICATIONS



Characteristics:




Cleanroom compatible



Pharmaceutical applications



Biohazard Applications



Stainless Steel



Autoclavable

For the following vacuum cleaner systems:

CWR-8

CWR-10

CWR-15

Please contact us for other systems.

- Compatible with ISO Class 4 (Former Class 10) Cleanroom conditions. Suitable for pharmaceutical applications.
- For the safe recovery and disposal of needles and glass shards.
- Light weight and easy to use, this unit will prevent needles and glass shards from cutting and getting stuck in suction hoses.
- This unit can be attached to models CWR-8, 10, 15 Cleanroom vacuum cleaner systems.
- Can be used with either the silicone Ø2" (50 mm) Tiger-Claw or Conical Nozzle attachments.
- All Stainless Steel AISI 304 construction.
- Completely autoclavable.
- 1.25" (32mm) wand connection.

Dimensions

Length:	21.5" (54.6cm)
Width	4 5/16" (11cm)
Height	3 3/8" (9.4cm)
Weight	3lbs (1.36kg)

Tiger-Vac®
...The name to vacuum with

Tiger-Vac Inc. (USA)
73 S.W. 12 Ave. Bldg. 1, Unit 7
Dania, FL 33004
USA
Tel.: (954) 925-3625
Fax.: (954) 925-3626
E-mail: sales@tiger-vac.com

Tiger-Vac International Inc.
2020 Dagenais Blvd. West
Laval, Quebec
Canada H7L 5W2
Tel.: (450) 622-0100
Fax.: (450) 625-3388
E-mail: support@tiger-vac.com

Toll free telephone: 1 800 668-4437 (Fax: 1 800 668-4439)

<http://www.tiger-vac.com>